



ArterioFlow® 7500

Arterial Compression Pump

Product Overview

The **ArterioFlow® 7500** is an arterial compression pump that provides powerful sequential compression to the foot and calf through sleeve garments to empty the veins in the legs and feet, permitting an increase in blood flow. Increased blood flow may help patients feel healthier and can allow them to move easier.

Features and Benefits:

- ✓ Convenient:** Adjustable therapy time, approved for in-home use.
- ✓ Powerful:** Effective compression for chronic disorders caused by reduced blood supply.
- ✓ Easy to Use:** LCD screen, simple setup.

Device Specifications:

Weight:	12 lbs.
Inflation:	Up to 4 seconds
Cycle Time:	Adjustable 15 - 30 seconds
Therapy Time:	Adjustable 5 - 180 minutes
Dimensions:	4" (H) x 10" (D) x 11.5" (W)
Electrical:	120 VAC, 60Hz, 100 Watts Max
Fuse Rated:	250VAC, 1.0AMP, SLO-BLO
Applied Part:	TYPE B
User Interface:	LCD Screen
Operation Mode:	Factory Default / Customer Set



Item Code: P-7500
LB05.0002 Rev. D 20160128