



CircuFlow™ 5200

Lymphedema Compression Therapy

Product Overview

The **CircuFlow™ 5200** is a compression therapy pump that creates intermittent gradient sequential pressure that is applied through the multi-chambered garments to increase venous blood flow and provide extracellular fluid clearances for the patient. It is used to help relieve painful swelling associated with lymphedema and other chronic venous disorders.

Features and Benefits:

- ✓ **Safe:** Device lock & easy shut off.
- ✓ **Custom Settings:** Two modes which can be set individually.
- ✓ **Flexible:** Able to skip chambers to prevent discomfort and customize cycle time.

Device Specifications:

Weight:	12 lbs.
Inflation Cycle:	Adjustable 5 - 20 seconds per chamber
Deflation Cycle:	12 seconds
Dimensions:	11.5" (W) x 10" (D) x 4" (H)
Cycle Time:	Adjustable 20 - 80 seconds
End Therapy Vacuum Time:	60 seconds
Therapy Time:	Adjustable 15 - 180 minutes
Electrical:	120VAC, 60Hz, 100VA MAX
Fuse Rated:	250VAC, 1.0AMP, SLO-BLO
Operation Mode:	Continuous Operation with Intermittent Loading



Item Code: P-5200
LB04.0002 Rev D 20160115