


Cirona[®] 6100

Deep Vein Thrombosis Prevention Therapy System

Product Overview

The **Cirona[®] 6100** Deep Vein Thrombosis (DVT) Prevention Therapy System is a pneumatic compression device that can help reduce the incidence of deep vein thrombosis, a potentially life-threatening condition in which a blood clot forms in a vein deep inside the body. The system consists of a pump and a pair of soft compression garments that wrap around the calf, calf-thigh, and foot on each leg or foot. Air compression inflates the garments to mimic the body's natural walking functionality, ensuring circulation, and prevent early blood clotting.

Features and Benefits:

- 
Reliable: Effective compression therapy, may reduce potentially fatal DVT.
- 
Customizable Two modes; single leg and dual leg for personalized treatment.
- 
Comfortable: Soft, disposable single use garments.

Device Specifications:

Weight:	5.5 lbs.
Dimensions:	4" (W) x 7" (L) x 11" (H)
Inflation Cycle:	12 seconds
Rest Period:	48 seconds
Cycle Time:	60 seconds
Pressure:	40mmHg
Electrical:	120 VAC, 60Hz, 45VA MAX
Fuse Rated:	250 VAC, 1.0AMP, SLO-BLO
Operation Mode:	Continuous



Item Code: 6100
LB37.0001 Rev G 20160115