




## Cirona<sup>®</sup> 6200

### Deep Vein Thrombosis Prevention Therapy System

## Product Overview

The **Cirona<sup>®</sup> 6200** is a compression therapy pump used for patients who may be at-risk of developing a painful and potentially life threatening blood clot known as deep vein thrombosis or DVT. The device comes with a pair of soft compression garments that wrap around the calf, calf-thigh, or foot on each leg or foot. Air compression inflates the garments during a 60-second automatically timed cycle (an approximate 12-second inflation period followed by a 48-second relaxation period) to mimic the body's natural walking functionality, to enhance circulation and prevent early blood clotting

## Features and Benefits:

- 
**Convenient:** comes with long-lasting rechargeable battery, bed hanger attaches easily in hospital or home care setting.
- 
**Effective:** Enhances circulation to provide therapy that may prevent clotting.
- 
**Comfortable:** Disposable compression garments made of soft, breathable material.

## Device Specifications:

<b>Weight:</b>	5.1 lbs.
<b>Dimensions:</b>	4" (W) x 7" (L) x 11" (H)
<b>Inflation Cycle:</b>	12 seconds
<b>Rest Period:</b>	48 seconds
<b>Cycle Time:</b>	60 seconds
<b>Pressure:</b>	40mmHg
<b>Electrical:</b>	100 - 240 VAC, 50X60Hz
<b>Fuse Rated:</b>	250 VAC, 1.0AMP, SLO-BLO
<b>Operation Mode:</b>	Continuous
<b>Battery:</b>	Lithium Battery, 10.8 V, 4400mAh (rechargeable)



Item Code: 6200  
LB37.0004 Rev C 20160301