

## Cirona<sup>®</sup> 6400

### Deep Vein Thrombosis Prevention Therapy System

## Product Overview

The **Cirona<sup>®</sup> 6400** is a portable Deep Vein Thrombosis Prevention therapy device. The device helps enhance proper circulation of blood in the body by mimicking the body's natural muscle contractions, thus increasing venous blood flow to the patient's extremities. The **Cirona<sup>®</sup> 6400** requires no set up and is one-button operation making it a convenient device for patients. It is fully portable and capable of delivering continuous DVT prevention for up to 24 hours for both acute care and home care settings.

## Features and Benefits:

- 
**Long Lasting Batteries:** Rechargeable batteries with up to **24 hours** of treatment
- 
**Portable:** Providing 100% mobility with no cords or tubing attached
- 
**Easy Operation:** Turn on/off with the press of a button
- 
**Adjustable Sleeves:** Soft, breathable fabric that adjusts to comfortably fit.

## Device Specifications:

<b>Weight:</b>	0.33kg / 0.723lbs
<b>Dimensions:</b>	70x 130 x 30 (mm) 2.75 x 5.1 x 1.2 (in)
<b>AC/DC Adapter:</b>	AC Input: 100~240 Volts, AC 50/60Hz, 10VA MAX DC Output: 5V 2A (2 Plugs)
<b>Operation Mode:</b>	Continuous / Up to 24 hours
<b>Battery Type:</b>	Lithium battery 3.7v 2250mAh (rechargeable)
<b>Charging Time:</b>	3 Hours
<b>Default Pressure:</b>	40 mmHg ± 5 mmHg
<b>Cycle Time:</b>	60 seconds



Item Code: 6400  
LB61.0003 Rev A 20170328